



**Wellington  
Hospitals  
Foundation**

## **Wellington Hospitals Foundation Favourite – Knitted Chevron Blanket**

FINISHED SIZE: 60 x 60cm in 4ply yarn



ABBREVIATIONS: md – mitred decrease (slip2, k1, pull slipped stitches over knitted stitch); ssk – slip, slip, knit (slip 2 stitches knitwise, then knit them together through the back loop); yo – yarn over.

WOOL: 4ply wool using 4 shades. You will need one ball of each, and 2 balls of the main colour.

### **INSTRUCTIONS:**

With first colour, cast on 193 sts.

Row 1: knit.

Row 2: k1, ssk, k9, yo, \*k1, yo, k10, md, k10, yo\* repeat\* to \*6 more times, k1, yo, k9, k2tog, k1.

Row 3: K, knitting into back of yo's.

Row 4: As row 2.

Row 5: As row 3.

Row 6, and every RS row: As row 2.

Row 7 and every WS row: Purl.

Keep repeating in colour sequence of your choice. We did 18 rows of each colour.

To finish off, do 4 rows of garter stitch as at the start.

**Thank you for knitting this blanket for Wellington babies**